


# OL 700-22

Pulsed Xenon Flash Source



The OL 700-22 Pulsed Xenon Flash Source is a high-pressure, short arc xenon flash lamp operating at 100Hz rate. Designed for applications requiring higher ultraviolet output, such as reflectance and transmittance and lower stray light than can be achieved with tungsten halogen source, it also offers lower spectral variation than tungsten halogen sources. The OL 700-22 comes with an SMA 905 connector that couples the source through a miniature fiber optic probe (optional) to accessories, including integrating spheres, transmittance, and reflectance attachments. It provides pulse-to-pulse stability, and the output stability improves with additional sample averaging. Pulse output stability can be improved by a factor of three by averaging as few as ten scans.

When the OL 700-22 is coupled to the OL 770 Multichannel Spectroradiometer, optically isolated signals are capable of synchronizing the operation of two light sources. Operation with a standard OL 700-22 source control cable permits single source operation. When the bifurcated OL 700-22 source control cable is used with two sources, each source can be independently controlled by the OL 770. The ability to cascade two OL 700-22 sources enables the user to perform near simultaneous reflectance and transmittance sample analysis.

Because it produces a pulsed signal, the OL 700-22 is less likely than continuous sources to contribute to solarization in optical fiber assemblies, which can occur when fibers are illuminated with signal wavelengths less than 260 nm.

#### Data Sheet: B090 Dec 2020 | Rev A

As part of our policy of continuous product improvement, we reserve the right to change specifications at any time.

### SPECIFICATIONS

<b>Spectral Range</b>		200 to 1100nm*
<b>Arc Voltage (V)</b>	Internally Preset	600 Volts
<b>Arc Current (I)</b>	$V \cdot (\text{Capacitance (F)} / L \text{ Inductance (H)})^{1/2}$	400 Amperes
<b>Discharge Energy (E)</b>	$.5 \cdot \text{Capacitance (F)} \cdot V^2$	.04 Joules
<b>Average Power</b>	Joules (E) * Pulse Rate (Hz)	4 Watts
<b>Pulse Rate</b>	Internally Preset	100 Hz
<b>Output Variation</b>	Day to Day (10 770UV Scans Averaged)	<2% (250 to 1100nm)
<b>Lamp Lifetime</b>	115 Days of Continuous Operation	1*10 <sup>9</sup> Pulses
<b>Timing</b>	Synchronized to OL 770 Exposure Signal	Single or Dual Operation
<b>Optical Interface</b>		SMA-905 Fiber Coupler
<b>Electrical Interface</b>	Opto-coupled	High Density -15 "D"
<b>Dimensions</b>		8" x 3" x 6" (20.3 cm X 7.6 cm X 15.2 cm)
<b>Weight</b>		2lbs. (900mg)

\*See Optional Fibers Section

### OPTIONAL FIBERS

<b>770-7G-RS</b>	VIS/ VIS/ NIR; 380 – 1100 nm; Glass FOP
<b>770-7Q-SMA</b>	UV/ VIS; 200 – 800 nm; Quartz Grade FOP

For more information visit [OptronicLabs.com](http://OptronicLabs.com)