

# TOTAL LUMINOUS FLUX (TLF) MEASUREMENTS OF LEDS AND OLEDs


## WHAT AND WHY

Total luminous flux (TLF) measurements are arguably one of the most common types of measurements used for characterizing light output from single die and multi-LED arrays, OLEDs and many other types of sources, such as incandescent, fluorescent, metal halide lamps, etc.

There are several different ways to present data from measured sources - radiometrically, photometrically, or spectroradiometrically. TLF is simply a photopically weighted measurement of the total energy emitted (*Watts*) from a device over the wavelength range 380 – 780 nm, referred to as the visible spectrum. The weighting factors are defined by the  $V(\lambda)$  function, which describes the response of the human eye in the visible region of the spectrum. Multiplying the radiometric (watts) measurement by the  $V(\lambda)$  function yields a measurement result in lumens (lm) and is most commonly referred to as TLF.

## HOW

There are two generally accepted methods of obtaining the total luminous flux of sources, goniometrically or using an integrating sphere. Each of these methods has their pros and cons.



## GONIOMETRIC MEASUREMENTS


- Are very time consuming.
- Requires data to be taken 360° in both the  $\theta$  and  $\phi$  directions.
- Equipment is costly, especially for large sources.
- A dedicated room is needed.
- Calibration can be tricky.

## INTEGRATING SPHERE

- Easy to use.
- Most widely accepted method for measuring TLF.
- Very fast, samples are placed at the center and data is acquired in a single scan.
- Simple to calibrate using OL 770's integrated calibration routine and aux lamp.

## SPECS/FEATURES:

- Convenient USB interface
- 25+ spectral scans/second
- Meets CIE 127 guidelines, Conditions A & B, TLF
- Auxiliary lamp for absorption correction and system throughput changes
- Low stray light performance
- High spectral resolution
- High sensitivity
- High dynamic range
- 0.5 nm wavelength accuracy
- Research-grade precision
- Compact, lightweight, portable enclosure
- Rugged strain relief and self-centering adapter



## IN SUMMARY:

- To make TSF measurements, either a sphere or a goniometer coupled to a spectroradiometer can be used, though the sphere is more efficient and is less susceptible to alignment and calibration errors.
- The DUT is placed inside the sphere or at the sphere wall (*i.e. the OL IS-1800 or OL IS-670-LED 6" Integrating Sphere*) and the OL 770 Spectroradiometer performs a system response calibration to compensate changes in throughput of the sphere due to the presence of DUT, holder, etc.
- The OL 770 Spectroradiometer records the Total Luminous Flux value and automatically calculates key parameters, such as Chromaticity, Lab/Luv, Peak Wavelength, Half Bandwidth, etc.

**OPTRONIC™**  
LABORATORIES

**Application Note: A20 Jan 2019**

As part of our policy of continuous product improvement, we reserve the right to change specifications at any time.