

NEW OL 83A vs OLD OL 83/83DS

	OL 83A	OL 83/83DS
CURRENT ACCURACY	± 0.01%	± 0.05%
LINE/LOAD REGULATION	± 0.01%	± 0.05%
COMPUTER INTERFACE	RS-232 (<i>Standard</i>) IEEE-488 (<i>Optional</i>)	None
FLEXIBILITY	Programmable	None
POWER REQUIREMENTS	115v/60hz or 230v/50hz	115v Only
PRICE	Less Than OL 83DS!	

New Features

- Simultaneous digital readout of lamp current, voltage & power.
- Ability to accept an operating target in terms of lamp current, voltage or power.
- Automatic calibration due warning.
- Lamp “library” for storing & recalling up to 10 user-defined lamp operating points with separate timers for each lamp.

Information Sheet: IS14 Jan 2019

As part of our policy of continuous product improvement, we reserve the right to change specifications at any time.