

OL 770 MULTI-CHANNEL SPECTRORADIOMETER RECOMMENDED COMPUTER CONFIGURATION

	OPTIMUM CONFIGURATION For optimum performance, use the following configuration:	MINIMUM CONFIGURATION The following configuration is adequate for OL 770 operation, but there will be noticeable decreases in response time and system performance:
COMPUTER	Intel® Pentium® III or Athlon®-Based IBM® -Compatible Personal Computer or Laptop, 1 GHz or Faster	Pentium® II-Based IBM®-Compatible Personal Computer or Laptop, 500 MHz
MEMORY	128 Mb RAM or Greater	64 Mb RAM
OPERATING SYSTEM	Windows® XP, Windows Vista™, and Windows 7	Same
MONITOR	SVGA Capable of Displaying 1024x768 or Greater in 24 Bit Color	SVGA Capable of Displaying 800x600 in 16 Bit Color
MASS STORAGE	Hard Drive > 1 Gb	> 20 Mb Free Space
MOUSE	Microsoft® Compatible	Same
PRINTER	Windows® Operating System Supported Printer	Same
INTERFACE (STANDARD)	USB – USB Port 1.1 or 2.0 /2	RS232 ~ 2 Spectra/Sec
SOFTWARE	Microsoft Vista™ and Windows 7	Same

* Use of the minimum system configuration will limit transfer rates in USB mode.