

SOFTWARE APPLICATION OPERATING SYSTEMS COMPATIBILITY CHART

PROJECT	OPERATING SYSTEM				
	WINDOWS XP	WINDOWS XP X64	WINDOWS 7	WINDOWS 7 X64	WINDOWS 8
OL 16A/AV, 65A/AV, 83A/AV SDK	Yes	No	No	No	Untested
OL 462 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 700-10 Software Development Kit	Yes	No	No	No	Untested
OL 700-30 Software Development Kit	Yes	Yes	Yes	Yes	Untested
OL 730C/CV SDK	Yes	Yes	Yes	Yes	Untested
OL 730D Software Development Kit	Yes	Yes	Yes	Yes	Untested
OL 730D Pulse Profiling Application	Yes	Yes	Yes	Yes	Untested
OL 750 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 750 Windows Demo Application	Yes	Yes	Yes	Yes	Untested
OL 750 Software Development Kit	Yes	Yes	Yes	Yes	Untested
OL 750-421S Software Development Kit	Yes	Yes	Yes	Yes	Untested
OL 756 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 756 Software Development Kit	Yes	Yes	Yes	No	Untested
OL 770 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 770 Software Development Kit	Yes	Yes	Yes	Yes	Untested
PCSRad Windows Application	Yes	Yes	Yes	Yes	Untested
PCSRad Software Development Kit	Yes	Yes	Yes	Yes	Untested
Keithley 24xx Software Development Kit	Yes	Yes	Yes	Yes	Untested
Lumesphere Plus Application	Yes	Yes	Yes	Yes	Untested
OL 490 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 490 Software Development Kit	Yes	Yes	Yes	Yes	Untested
HSi U-MX Windows Application	Yes	Yes	Yes	Yes	Untested
OL 731 Windows Application	Yes	Yes	Yes	Yes	Untested
OL 700-MS Application (<i>Modular Software</i>)	Yes	No	Yes	No	Untested

Information Sheet: IS29 Jan 2019

As part of our policy of continuous product improvement, we reserve the right to change specifications at any time.