

SPECTRORADIOMETER IMPROVES ACCURACY OF PHOTOVOLTAIC MEASUREMENT

Reprinted with Permission from Laser Focus World,
PennWell Corp.

January 2008

Spectroradiometer improves accuracy of photovoltaic measurement

GHASSAN E. JABBOUR, EVAN WILLIAMS,
RICHARD YOUNG, AND ALEXANDRE Y. FONG




FIGURE 1. The current response for a commercially available amorphous-silicon photodiode as a function of illumination follows an exponential curve.

S

ince the discovery of the photoelectric effect in 1939 by Alexandre Edmond Becquerel, the performance of photovoltaic (PV) detectors has gradually improved from scientific novelty to a potentially viable alternative energy source. Photovoltaics are semiconductor devices that convert light into electrical current through the creation of electron-hole pairs at a *p-n* junction forward-biased by incident photons with energy greater than its bandgap energy.¹ Conventional silica-based devices—which comprise the majority of available devices—have conversion efficiencies of less than 15%. Other materials such as cadmium telluride have shown laboratory efficiencies of almost 20% and even 30%, but are limited in practical application.

Until recently, photovoltaic solar energy

has been too costly. However, precipitated by international events including increased global demand for energy and


concerns about limitation of petroleum sources, recent interest in alternative and renewable energy sources has refocused attention on and investment in the potential of PVs. Active programs are under way at the U.S. Department of Energy's (DOE) National Center for Photovoltaics at the National Renewable Energy Laboratories (NREL; Golden, CO) and by large industrial energy firms such as British Petroleum's (London, England) Solar Products group.

While it is impossible to reproduce the brightness level of the Sun, it is possible to simulate a photovoltaic detector's response with improved accuracy by integrating a white-bias source into a scanning monochromator-based spectroradiometer.

BIOGHASSAN E. JABBOUR is a professor of the School of Materials and director of the Center for Advanced Photovoltaics, and EVAN WILLIAMS is a graduate student in the Optoelectronic Materials and Devices group at the Flexible Display Center, both at Arizona State University, 1290 S. Normal Ave., Tempe, AZ 85284; e-mail: jabbour@asu.edu. RICHARD YOUNG is vice president of science and technology and ALEXANDRE Y. FONG is vice president of sales and marketing at Optronic Laboratories, 4632 36th St., Orlando, FL 32811; e-mail: afong@olinet.com; www.olinet.com.

Further evidence of this renewed interest is the entry of other industry leaders into the fray such as Applied Materials (Santa Clara, CA; see www.laserfocusworld.com/articles/274723) and record attendance at the recent European Photovoltaic Solar Energy Conference and Exhibition in Milan, Italy.

According to the DOE, the worldwide market for PVs reached \$4.7 billion in 2003, and according to the 2005



Polarizers

UV - visible - near IR

- High transmission
- Excellent extinction ratio
- Custom shapes and sizes
- Experienced OEM supplier
- Mounted and unmounted
- Standard materials/
custom designs
- Laser and imaging
instrumentation applications

Call for a free copy of our new
2007-08 Polarization Handbook.

Telephone: 303.833.4333
E-mail: mlo@meadowlark.com
Website: www.meadowlark.com

meadowlark optics
polarization solutions

PHOTOVOLTAICS, continued




FIGURE 2. The current reference spectral irradiance for photovoltaic (PV) testing is the ASTM G173-03 reference spectra. (Courtesy of ASU)

Solar Photovoltaic Report from ABS Energy Research (London, England), the PV market is expected to triple over the next three years. Significant contributions to this growth will come once the fabrication of efficient and reliable PVs becomes a commercial reality.²

New and composite semiconductor materials, polymers, thin-film coatings, and nanotechnology such as quantum dots all show tremendous promise toward the goals of improved conversion efficiencies, lower installation costs, and adaptability to broader applications for PV solar-cell detectors. Critical to the work of improving PV conversion efficiencies is determining and expanding their performance across a wider spectral range. For example, developments in new polymers that absorb light from the ultraviolet through the infrared offer unforeseen and exciting possibilities in terms of conversion efficiencies and applications.³ Accurate characterization of PV spectral response is central to this objective.

Solar-cell performance and power efficiency

The primary parameter describing the performance of a solar cell is the power-conversion efficiency, as evidenced by the watt-peak rating given to modules. The most standard test used to characterize a solar cell is the current-voltage (I-V) response under illumination, and from this, the power-conversion efficiency of the device can be calculated (see Fig. 1).

There are several other PV parameters of importance. The short-circuit current density, J_{sc} , is a measure of the photocurrent with no external bias applied and the two electrodes connected via an external circuit. The open-circuit voltage, V_{oc} , rep-

resents the necessary applied voltage for which no net current passes through the device. There is a point along the response curve for which the product of the voltage and current gives a maximum; this is the point of maximum power and represents the load the cell should drive for optimal output. A rectangle can be created using the origin, the point of maximum power and perpendicular lines to the

axes (shaded box). The ratio of areas for this rectangle and one formed using the origin, V_{oc} , and J_{sc} points (dotted box) is known as the fill factor and is calculated as $P_{max}/J_{sc}V_{oc}$. The power conversion efficiency η can be expressed as:

$$\eta = V_{oc} \cdot J_{sc} \cdot FF / P_{inc}$$

Here, P_{inc} is the incident optical power. It is therefore of great importance to accurately know the incident power to calculate the correct efficiency.

Furthermore, in addition to the illumination source intensity, it is critical to account for the spectrum of the illumination source and its difference from the standard solar spectrum to obtain the "true" power efficiency. A reference cell is most commonly used to establish the intensity of the illumination source used for testing. There may be an error in the measured J_{sc} of the test device arising from differences in the spectrum of the illuminator and/or differences in the spectral response of the reference cell and the device under test. A spectral mismatch factor can be calculated to express this error and is given by:

$$M = \frac{\int_0^\infty SR_R(\lambda) E_{STND}(\lambda) d\lambda}{\int_0^\infty SR_R(\lambda) E_s(\lambda) d\lambda} \cdot \frac{\int_0^\infty SR_T(\lambda) E_s(\lambda) d\lambda}{\int_0^\infty SR_T(\lambda) E_{STND}(\lambda) d\lambda}$$

Here, SR_R is the spectral response of the reference cell, E_{STND} is the reference spectral irradiance, SR_T is the spectral response of the test cell, and E_s is the test-source spectral irradiance.

The current spectral irradiance reference is the ASTM G173-03 reference spectra (see Fig. 2). Using a properly characterized reference cell and knowing the test-source spectral irradiance and test-device spectral response allows

for accurate determination of the power efficiency for a solar cell. Considerable efforts have been made by the international community to create standard testing methods to provide fair analysis of device performance from different researchers.⁴

Basic spectral-measurement methodology

Normally, detector spectral response is measured by exposing the detector to monochromatic light that is scanned through the wavelength range of interest. The amount of output signal resulting from exposure to the monochromatic light at each wavelength increment is first measured with a standard detector and then the test detector is placed in the light. The signal from the test detector divided by the amount of light is the response of the detector. The source is typically a scanning grating monochromator (see Fig. 3).

A scanning monochromator uses precision mechanics and slits to select specific wavelengths dispersed by gratings from a broadband spectrum. However, the resulting amount of monochromatic light from such a source is very small. It

has also been observed that many solar cells respond differently at high light levels and low light levels. While it is impossible to get the monochromatic light the same brightness as the sun, a white bias light source can be used to effectively simulate a bright environmental condition. This is just a constant light source that raises the overall light level to the same value as the sun (see also www.laserfocusworld.com/articles/286515). In this way, a more representative measure of the PV solar cell's spectral performance can be made.

White-light biasing of solar cells

The bias source is provided to illuminate the test device with approximately 1 solar constant (1370 W/m^2) to simulate the actual operating environment. Utilization of a bias source will require an AC




FIGURE 3. A scanning grating monochromator is typically employed as the source when performing detector spectral response measurements. The internal layout of a double monochromator shows the broadband source, collimated output, and reference detector and detector under test. (Courtesy Optronic Laboratories)

signal lock-in detection system. The bias source is shuttered when performing the system calibration with the reference detector and then used to bias the test detector and PV solar cell during the measurement. The test detector will respond to this DC light source and the chopped monochromatic beam. An additional


Green YAG Lasers for Solar Cell Manufacture

Since 1996 our 532-nm Nd:YAG laser sources have been used to manufacture PV products.

Lee Laser Nd:YAG lasers at 532 nm set the standard for:


- High-power lasers to process polycrystalline silicon
- Low-power lasers for fine, thin-film structuring
- Up to 300 Watts
- Superior beam quality
- Industrial grade reliability

**LEE
LASER**




leelaser.com

Orlando, Florida USA Phone: 407-812-4611
Email: salesdept@leelaser.com




OVER
2,261
NEW OPTICAL PRODUCTS

TECHSPEC®




MECHANICS



BeamX™ BEAM
EXPANDERS



F THETA LENSES



LASER MIRRORS

MORE CHOICES
FASTEST DELIVERY

NEW!
404-page
Master Source
Book

Request your
FREE copy today!

800.363.1992

www.edmundoptics.com/LF

ED Edmund
optics | worldwide

Visit us at Photonics West, Booth 1019

PHOTOVOLTAICS, continued




FIGURE 4. The relative response of a PV solar cell at 980 nm without the white-light bias (top) is contrasted to the response of the same cell at 980 nm with a 1 kW/m² white-light bias (bottom). (Courtesy Optronic Laboratories)

detector preamplifier (OL 750-613 from Optronic Laboratories) decouples the DC response of the detector to the bias source so that only the AC component of the response is measured.

Spectral-response measurements

An Optronic Laboratories OL 750D double monochromator spectroradiometer with a white-light bias was used to measure the spectral response of PV cells. Since the bias light is DC and the monochromatic light is AC chopped, the detector response can be measured by locking in to only the AC portion of the detector signal. The white-light bias is a small spot, typically 1 mm in diameter, to ensure that the detector or amplifier is not saturated. The monochromator light is focused to a spot about half the size of the white-light spot and completely within it. The white-light bias level can be set by adjusting the current to the white light lamp.

It is important that the white-light bias coincides with the monochromatic light at the sample. If they do not coincide, the same responsivity is seen in both biased

and unbiased conditions. Because the spots are small, alignment is critical. Also, many polycrystalline materials exhibit nonuniform responsivity over the wafer. A small spot is ideal for characterizing the uniformity of such devices (see Fig. 4).

Reflectance and efficiency measurements

Determining the internal quantum efficiency of a PV solar cell requires measurements of the reflectance of the device. Reflectance from such cells is typically composed mainly of specular and some limited diffuse component. Although the diffuse reflectance can be small when compared to the specular component, it can have a large effect on the device quantum efficiency. Adding an assembly to measure the specular and diffuse reflection components such as the OL 750-71 provides measurements of both power response and reflectivity of the test detector, allowing the calculation of external and internal quantum efficiency. The reflectance measurement geometry is suited to specular (mirror-like) samples; the OL 740-70 diffuse reflectance accessory should be used for samples without mirror-like finishes.

Despite continuing polysilicon supply issues, the PV industry continues to grow and mature, with a variety of different materials in the works to further improve conversion efficiencies—in fact, plastic solar cells have been reported with efficiency values exceeding 6% (see www.laserfocusworld.com/articles/294199). All of this underscores the need for accepted, accurate, and meaningful metrology of solar-cell spectral performance evaluated in realistic operating environments. White-light biasing and high-performance monochromatic sources are fundamental to this objective. □

REFERENCES

1. M. A. Green, *Energy Policy* 28, 989 (2000).
2. S. E. Shaheen et al., *Appl. Phys. Lett.* 79, 2996 (2001).
3. R. Janssen, *SPIE Newsroom* (online), 10.1117/2.1200606.315 (2006).
4. K. Emery and C. Osterwald, *Solar Cells* 17, 253 (1986).

Tell us what you think about this article. Send an e-mail to LFWFeedback@pennwell.com.